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
Fig. on front page: "Main" display of a DCU system for FlexAct® COM

Sartorius Stedim Biotech GmbH

Technical Documentation | Technische Dokumentation
August-Spindler-Strasse 11

D-37079 Goettingen

Germany | Deutschland
Internet: www.sartorius-stedim.com
E-mail: tech.pubs@sartorius-stedim.com

Printed in Germany on paper that has been bleached without any use of chlorine.
If you require mandatory information regarding specific properties of your
FlexAct® COM that goes beyond that which is described here, please contact us.

No part of this publication may be reprinted, reproduced or translated in any form or
by any means without the prior written permission of Sartorius Stedim Biotech GmbH.
FlexAct® is a registered trademark of Sartorius Stedim Biotech GmbH

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Contents

Introduction	5	7. Phases, Step Controls	35
1. Start-up of the DCU System	9	7.1 Overview on Phases	35
1.1 General Notes	9	7.2 Phases of the FlexAct® System	35
1.2 Functional Groups and Submenus	9	8. Main Function "Settings"	39
1.3 System Behaviour at Start-up	10	8.1 Overview	39
2. Basics of Operation	11	8.2 Setting of System Parameters	41
2.1 Touch Display "Main"	11	8.3 Setting of Process Value Ranges	42
2.2 Appearance of Touch Keys	13	8.4 Manual Operation	43
2.3 Overview on the Main Function Keys	14	8.5 Menu "Service"	54
2.4 Direct Function Keys for Access to Submenus	15	9. Supplement	55
2.5 Lists and Tables	16	9.1 Alarms and Messages	55
3. Main Function "Main"	17	9.2 Trouble Shooting for the DCU system	57
3.1 Features and Functions	17	9.3 Locking Functions	57
3.2 Functional Elements of the "Main" Display	18	9.4 Abbreviations	58
4. Main Function "Trend"-Display	19	9.5 Arrangement and Assignment Plans	59
4.1 Intended Use	19	9.6 Password System	61
4.2 Menu Overview of the Main Function "Trend"	20		
4.3 Setting of the "Trend" Display	21		
5. Main Function "Calibration"	23		
5.1 Main Display "Calibration"	23		
5.2 Operation of the Calibration Function in General	24		
5.3 pH Calibration (Conventional pH sensors)	25		
6. Main function "Control Loops"	29		
6.1 Overview on Controllers	29		
6.2 Controller Operation in General	30		
6.3 Parameterization of Controllers in General	31		
6.4 Pump Controller	33		

This page intentionally left blank.

Introduction

Intended Use, Design and Functions of the Device

The DCU system is a digital measuring and control system for controlled devices, like bioreactors, crossflow units and other laboratory equipment, e.g. FlexAct® systems or incubation shakers, of Sartorius Stedim Biotech GmbH. The DCU functions described herein bases upon the DCU4 system for bioreactors and are especially adapted to the FlexAct® systems.

The DCU hardware can either be contained in the controlled device or in a stand-alone table top housing or be designed as a rack version in control cabinets of pilot and production plant bioreactors.

The measurement and control functions are implemented in the software. Standard functions include the measuring of parameters important for the intended process, the monitoring of their run, the calibration of sensors and dosing counters and corresponding controls associated with the process parameters. The extent of functions depends on the controlled device and can be adapted to the needs of the customer.

DCU systems can be integrated into hierarchical automation systems. A host computer can be connected, which is running the proven MFCS-System, for instance, providing functions like graphical visualization of processes, data acquisition and storage, process protocols, etc..

The digital control system DCU is an example of the product program of Sartorius Stedim Systems GmbH of bioreactors and peripheral equipment. Further information about the product program will be available from your representative of Sartorius Stedim Biotech or directly from:

Sartorius Stedim Systems GmbH
Robert-Bosch-Straße 5 – 7
D-34302 Guxhagen, Fed. Republic of Germany

phone +49.5665.407.0, fax +49.5665.407.2200
e-mail: info.fermentation@sartorius-stedim.com
WebSite <http://www.sartorius-stedim.com>

Documentation Release Notes

This is an "original" Operating Manual. It refers to standard functions of the DCU system and their adaptation to the FlexAct® systems, as far as known at the release date.

Displays, parameters and settings shown in this Operating Manual are only exemplary. Unless especially indicated they do not refer to and are not intended for the configuration and operation of a DCU system assigned to a specific controlled device.

This document, however, can describe functions, which are not included in a configuration delivered, or the functions can be implemented in a different way. Customized modifications can be described in supplemental documentation, which is added to this manual or supplied extra.

The controlled devices, their setup, equipment and operation during a process are described in additional documents. Information on the hardware of the DCU system and the specifications of the configuration delivered are included in the "Technical Documentation". Complementary documents can be part of the "Technical Documentation" folder delivered with the device or be available extra.

Legend

Information provided herein which is of special importance and | or need be observed and followed with special care, is highlighted by its formatting. In particular, the highlighting follows these rules:

Warning notes indicate hazards, which can result in serious injuries to persons or death, at estimated medium risk.
Warning notes bear this symbol and are highlighted as shown here.

Caution notes indicate hazards, which can result in possible injuries to persons, at estimated minor risk.
Caution notes bear this symbol and are highlighted as shown here.

This symbol indicates electrical hazards (by voltage or current). Information on their cause and preventive measures are highlighted as shown here.

Notices indicate hazards, which will or may lead to damages to equipment, if the safety information is ignored.

This symbol indicates important aspects which need be considered, or operations which shall be executed with special care, in order to avoid erroneous handling or malfunctions.

▷ The "button" indicates other noticeable information.

1., 2., ... Numbered paragraphs show steps of operation in recommended order.

► The arrow in brackets refers to information which can be found elsewhere.

Safety Information

The safety and hazard notes given in this document are meant exclusively for the equipment, as delivered. They will complement other regulations related to the application, which the operators must observe in place. Applications which entail specific hazards can require additional safety equipment or measures for the protection of staff and work environment. This document neither considers such circumstances nor does it deal with specific safety regulations, such as legal or elsewhere compulsory rules.

Possible electrical hazards. Electrical shocks are possible at

- mains supply of wrong specifications or without proper grounding
- splashing of water or medium against energized mains sockets
- opening the cabinet while connected to the energized power supply.

Check the type label for the correct power supply. Ensure proper grounding of all equipment and splash-proof installation of the mains. Do not use multi-outlet power strips placed on the laboratory table. Keep all units switched-off during connection.

Opening of the control cabinet is not required for normal operation; it is restricted to authorized service staff for service works, only.

Leakage or damages of bags, tubings and connectors will release the medium, which can be hazardous (e.g. toxic or pathogenicous).

Provide a working place of appropriate safety class and equipment.

Risk of personal injuries.

Depending on the equipment, FlexAct® units can be very heavy.
Avoid moving over inclined floor or ramps.
Place at even surfaces. Carefully locate the stoppers.

Measured values and settings shown herein are default or exemplary.
Only if indicated they are necessary or suggested for operation.
Admissible settings of the controlled device and, if applied, customized specifications are provided in the configuration documentation.

Persons responsible for the application must check if hazards can occur during the process, either at permitted settings or erroneous operation. Users must be instructed carefully.

Only system administrators or authorized, well trained and experienced operators are allowed to change the system configuration.

This page intentionally left blank.

1. Start-up of the DCU System

1.1 General Notes

1. Start-up of the DCU System

After switching-on the control unit the DCU system boots and initializes the control system. This will take about 1–2 minutes. The system is ready for operation when the "Main" menu appears on the screen.

The "Main" menu provides a schematic overview of the controlled device. It is the starting point for operation:


Fig. 1-1: "Main" display of the Control Unit of a FlexAct® COM system

1.2 Functional Groups and Submenus

The operating functions of the system are organized as "Main functions":

- ▷ **"Main"**: Schematic overview on the controlled system with display of measured and calculated process values. Touch keys in the tags of the parameters provide access to commonly operated functions.
In the case of "Twin" versions the system starts with the overview "All" which contains the most important parameters and related submenus. The displays "Main 1 | Main 2" show their full extent.
- ▷ **"Trend"**: Function for visual monitoring of parameters.
- ▷ **"Calibration"**: Submenus for operating the calibration routines.
- ▷ **"Settings"**: Factory settings, manual setting of measured values and control signals for simulation, reset of the control system, etc.
- ▷ **"Remote"**: Enables | disables remote operation.
- ▷ **"Alarms"**: Access to overview display of alarm messages.

The operating information provided herein refer to standard functions. Details on the specific setup of a controlled device are provided with the configuration documentation or, if applied, the Functional Specification.

Prior to attending the DCU system, e.g. for data input, setting of functions, etc. note the information on the system behavior below.

1.3 System Behaviour at Start-up

After switching on and program start, as well as at return of the mains after a power failure of the control unit, the measurement and control system starts (restarts) with a defined basic state [→ configuration]:

- ▷ The configuration will be loaded:
 - initially and after a reset of the system as preset in the factory
 - during use with the parameters as stored in the meantime.
- ▷ Parameters set by the operator are stored in a battery buffered memory, they can be applied for the (next) process again:
 - Setpoints
 - Calibration parameters
 - Phases (if implemented in the configuration).
- ▷ Controllers are "off", actuators (pumps, valves) in default position.
- ▷ Timer sequences (e.g. phase controls) are "off" (stopped).

The start-up behavior of outputs and system functions which directly act on the controlled device (e.g. controllers, timer etc.), depend on the cause and duration of the break. These kinds of interruption are distinguished:

1. Switching "off" and "on" via the mains switch of the control unit or controlled device.
2. Shut down due to mains failure of the laboratory mains.

The operator can set a maximum power down time ("Failtime") in the submenu "System Parameters" of the "Settings" menu:


Fig. 1-2: Submenu "System Parameter", see description in section "Settings"

In the case of a power failure which is shorter than "Failtime", the system continues as follows:

- ▷ A message "Power Failure" displays the failure time and duration.
- ▷ Controllers run on with the setpoint active at the power failure.
- ▷ Timers and setpoint profiles are continued (if implemented).

If a power failure takes longer than "Failtime", then the control system behaves as if being switched-off, i.e. the system is set to its basic state.

After restart the system displays the message "Pwf stop ..." [► "Alarm messages" in the supplement], with date and time of the power failure.

2. Basics of Operation

2.1 Touch Display "Main"

2. Basics of Operation

The display "Main" provides a graphical overview on the setup of the controlled device with symbols of equipment included and peripherals connected, such as the FlexAct® module with DCU system, pump, filter modules, LevMixer® Palletank®, receiving bags (Palletanks®) and test points.

Headline with display of system state →

Working area:

- ▷ Display of functional elements*):
 - pH, Temp or Weight: measuring devices, e.g. showing measured values; access to menus for setting the alarm monitoring;
 - PUMP: direct access to menu for entering the setpoint and operating mode and submenu for setting and parameterization of the control function;
 - Filtration: access to submenu for running the filtration phase.

Bottom section with main function keys →

- ▷ Access to main displays of main functions


Fig. 2-1: "Menu Main" of DCU for FlexAct®

*) Schematic setup, functional elements, TAGs and parameters depend on the configuration.

2.1.1 Headline

The headline serves only for display of status information:

Fig. 2-2: Display of active main function


Selected Main function: example "Main"

2010-01-18 09:24:49

Date in format [yyyy-mm-dd]; time in format [hh:mm:ss]

Access to submenu for entering the process time (process start reference)
- available in selected configurations (e.g. if ordered)

Display of alarms:

- white bell: no alarm
- red bell: alarm occurred, details are provided in the alarm TAG
[▶ list of alarm tags in the supplement] and in the menu "Alarm"
- Alarm messages are available in the Main function "Alarm"

2.1.2 Working Area

- ▷ The "working area" displays the functional elements and submenus of the active main function, e.g. in case of the FlexAct® system:
 - Function keys with current process or set values
 - Pumps or scales with measured values
 - Controls, e.g. for "Pump", with current set value
 - Sensors, e.g. for "Temp", "pH", with their measured values
 - Peripheral equipment, e.g. LevMixer® Pallettank® and receiving bags (Pallettanks®).


Fig. 2-3: Example: Working area of "Main" display

- ▷ In general the DCU system is operated directly at the display. The functional elements in the working area and the "Main Function" keys in the bottom section are designed as "touch keys". Tapping the touch keys opens related submenus, e.g. as required for input of data or setpoints, or for selection of operation modes.
- ▷ Displays and available functions, "Tag" names, the parameters and the submenus depend on the controlled device to which the DCU system is assigned to, and on the implemented configuration.

2.1.3 Bottom Line

The "bottom line" has "Main function keys" for switching over between the main functions:

Fig. 2-4: Selection of the main function "Main" via its main function key

Appearance:

- Main function selected: button light-grey, sunk-in
- Unselected functions: buttons dark grey, elevated

2.2 Appearance of Touch Keys

The appearance of the touch keys, i.e. how the functional elements are displayed, indicates their current state and intended use:

Symbol	Display	Functional element, meaning, use
	Functional element	Field for short name („TAG“) of process value or controlled value, e.g. TEMP, pH, WEIGH, PUMP, Filtration /
	← touch key underlined with grey frame	Field for measured or controlled value in its physical unit: – submenu or function selectable by tapping the key
	← touch key underlined with green frame	Measuring or output of the functional element is active, the measured value or control variable is displayed: – frame light-green if controller is in cascade mode.
	← touch key underlined with yellow frame	Key display if function is switched to "man" mode: – can be either "manual off" or "manual on"; – no automatic control possible
	no frame below	For display only, no submenu assigned (function cannot be selected)

Fig. 2-5: Appearance of the Touch keys

- For examples on functional elements, their TAGs, measured values, parameters and submenus which can be selected via the "Touch keys" see the [► section "Main-Displays"].

2.3 Overview on the Main Function Keys




Key, symbol	Meaning, use
Main function "Main" 	Start-up display providing a graphical overview on the controlled device: <ul style="list-style-type: none"> – display of the components of the current configuration – overview on measured values and process parameters – direct access to submenus most important for operation.
Main function "Trend" 	Monitoring the course of a selection of 8 parameters from: <ul style="list-style-type: none"> – process values – setpoints of control loops – outputs of controllers (control loops).
Main function "Calibration" 	Menus for calibration functions. e.g. of: <ul style="list-style-type: none"> – Sensors for measuring pH
Main function "Controller" 	Operation and parameterization menus of controllers: <ul style="list-style-type: none"> – not available for current configurations DCU system for FlexAct® – pump control submenus accessible from "Main" display.
Main function "Phases" 	Main function on programmed routines (batches, timer, step controls): <ul style="list-style-type: none"> – not available for current configurations DCU system for FlexAct® – operation of filtrate phase accessible from "Main" display
Main function "Settings" 	Basic system settings, such as: <ul style="list-style-type: none"> – Measuring ranges of process values – Manual operation, e.g. of inputs, outputs, controllers, etc. – Communication with remote units (e.g. printers, host computers) – Configuration management (by authorized service only).
"Remote" operation 	Operation by remote (Host) computer systems: <ul style="list-style-type: none"> ▷ Tapping the "Touch-key" just enables the remote control; refer to section "Settings" for the configuration of the communication.
"Alarm" display 	Overview table of alarms which occurred during the process: <ul style="list-style-type: none"> ▷ At alarms the bell changes its color. ▷ Symbol displayed red : the table has unacknowledged alarms. ▷ Tapping the key opens the alarm table.

Fig. 2-6: Overview and appearance of the Main Function keys

All main functions can be selected at any time during a running process. The title of a selected main function, which is displayed in the working area also appears in the headline.

2.4 Direct Function Keys for Access to Submenus

- ▷ The functional elements displayed in the menu "Main" can include direct function keys providing direct access to the most important submenus or subfunctions for:
 - numerical inputs of setpoints, delivery rates, alarm limits, etc.;
 - selection of operation modes of controllers;
 - settings of or proceeding through calibration functions.

Which submenus or functions can be accessed directly from the main menu depends on the configuration. You can tap on the touch keys at the display to find out about your setup.

- ▷ The section "Main"-Display below has examples of the submenus (subfunctions) which can be accessed via the direct function keys. More detailed information on the associated functions and possible inputs is provided in the sections "Calibration" resp. "Controllers".

Example: Input of the setpoint of the transfer pump:

1. Tap the touch button of the "PUMP" symbol in the working area of the "main" menu.
- ▷ A submenu pops up, which has a numerical key pad for data entry on left hand and with buttons for "Mode" selection on right hand.
2. Enter the new setpoint (observe the permissible range within "min." and "max."). Delete wrong inputs and quit the submenu with the "C" button or reenter the setpoint again.
3. Confirm with [ok]. The submenu is closed. The new setpoint will be applied and displayed.


Fig. 2-7: Direct access for inputs and selection of the mode of the "TEMP" controller.

Example: Input of operation mode for the controller:


1. Tap the touch button of the "PUMP" symbol in the working area of the "main" menu.
2. Tap the touch button of the intended controller mode in the pad on right hand.
3. Confirm with [ok]. The submenu is closed. The new mode will be applied and displayed.

Additional Notes

- ▷ You can move forward via the key  to the detailed display of the PUMP controller.

2.5 Lists and Tables

If submenus contain lists of functional elements, tags or parameters which cannot be displayed in one window, the window has a slide bar on right hand with a "slider" key.

Ch.	Proc. Value	Min	Max	Alarm	Alarm Low	Alarm High
1	PRESS	0.000 bar	4.000 bar	enabled	0.000 bar	1.000 bar
2	WEIGHT	0.0 kg	1500.0 kg	disabled	0.0 kg	1500.0 kg
3	TEMP	0.0 °C	150.0 °C	disabled	0.0 °C	150.0 °C
4	pH	2.00 pH	12.00 pH	disabled	2.00 pH	12.00 pH
5	COND	0.0002 mS	1000.0000 mS	disabled	0.0002 mS	1000.0000 mS
6	PUMP	0.0 %	100.0 %	disabled	0.0 %	100.0 %
7	FLT_WGHT	0.0 kg	1500.0 kg	disabled	0.0 kg	1500.0 kg
8	SET_WGHT	0.0 kg	1500.0 kg	disabled	0.0 kg	1500.0 kg

Fig. 2-8: Selection of process values in the main function "Settings"

To move through lists or tables which include more entries than a single window can display:

1. Tap the arrow key down (▽) resp. up (△).
2. Tap and move the slider (light grey section in the slide bar) itself.
3. Tap the line at the slider bar at the relative position of the list or table, where the TAG of the functional element is arranged.

3. Main Function "Main"

3.1 Features and Functions



Fig. 3-1: Access to submenus from main menu via direct function keys (elements, tags and parameters depend on the configuration).

3.2 Functional Elements
of the "Main" Display

Functional element	Display	Touch-key function
	Parameter key showing the current mode or state: <ul style="list-style-type: none">– grey line: "auto - off"– green line: "auto - on"	Access to submenus for <ul style="list-style-type: none">– setting the alarm limits of measured PV's<ul style="list-style-type: none">– Temp– pH– Weight– setting of controls, e.g.:<ul style="list-style-type: none">– pump setpoint and operating mode– running the pump in reverse mode– running the filtration phase

*) Implemented functional elements, TAGs and parameters depend on the configuration delivered.

The overview on the previous page provides an example, which submenus are accessible by enabling the functional elements on the display "Main".

4. Main Function "Trend"-Display

4.1 Intended Use

4. Main Function "Trend"-Display

The "Trend" display provides a graphical display of the course of selectable parameters, such as process values, over a certain period of time. The resulting overview allows for a quick evaluation whether the process runs as intended or if, when and how it is affected by particular operator settings of parameters or by control actions or events. The period of time starts at enabling of the DCU system. The "Trend" display covers the run of the selected parameters for the preset protocol time.

Trends are not stored. For recording the run of parameters the data must be transmitted to a remote printer or host device (e.g. a SCADA system, such as the MFCS)!

- Selectable are parameters from either process values PV, setpoints or controller output signals, as implemented in the configuration:

For changing existing settings you can select "Channel # Settings" – "Select Buffered Channels"; with "#" as the no. of the channel [➤ example on next page for selection of the channels and viewing the parameter table].

- Possible display spans are the full range of each parameter and free selectable increments of the full range in terms of its physical unit [➤ "Channel # Settings" – "Min/max displayed value"].
- Available protocol times are 1, 12, 24, 36 or 72 h [➤ submenu "Display Range"].


Fig. 4-1: Trend display for monitoring the course of selected parameters over a given period of time (example without record of data).

4.2 Menu Overview of the Main Function "Trend"



Fig. 4-2: Overview on possible settings of the main function "Trend"

- ▷ Selectable parameters depend on the configuration. The setting shown for Channel 1 "TEMP" is exemplary for all other channels.

4.3 Setting of the "Trend" Display

Selecting parameters for the Trend display:


1. Enable the main function "Trend".
2. In the headspace tap the channel key "#" whose parameter should be set or changed. The submenu "Channel # Settings" appears.
3. To select another parameter for this channel, tap the touch key "PV".
 - ▷ The menu "Select Buffered Channel" displays the currently assigned parameters.
4. Tap "Cfg" to display the parameters of the current configuration.
5. If the intended parameter is not visible, you can move through the table with the arrow keys or the slider [▷ fig. "Overview", above].
6. Tap the touch key of the parameter and confirm the selection with [ok].
 - ▷ To unselect a parameter from the display at all, you can tap the touch key [.....].

Changing the display span of a parameter:

1. Select the window "Channel # Settings" and tap the key "Max" or "Min" which should be set.
2. Enter the upper (or lower) display limit within the given span. The full span is shown below the window for data entry.
3. Confirm all settings with [ok]. If the changes should be ignored, quit the window with "C".
4. Close the window "Channel # Settings" by tapping [X].

Reset of the display range:

- ▷ The display range of a selected parameter can easily be reset by tapping "Reset Range" in the window "Channel # Settings". This resets a modified setting of "Min" and "Max" to the default display range of the parameter.

Setting the color of a displayed parameter:

▷ The color is selectable from a set (table) of given colors.

1. Select "Channel # Settings" of the parameter which should have different color and tap the touch key with the name of the color.
2. In the table of colors tap the key with the name of the color, which should be applied. The selection will be assigned and enabled immediately.

Selecting a new „Display Range“:

1. Tap the touch key "

5. Main Function "Calibration"

5.1 Main Display "Calibration"

5. Main Function "Calibration"

The main function "Calibration" combines all calibration functions of configuration of a DCU-System, e.g. the calibration of sensors. Tapping the main function key "Calibration" opens a schematic overview of the controlled device, similar to the "Main" display, in this case with all functional elements relevant for calibration. The submenus of the calibration routines can then be accessed simply by tapping the assigned touch key.

Which calibration functions are available depends on the configuration. In case of the FlexAct®-system this is, e.g.:

▷ "pH": calibration of pH-sensor


Fig. 5-1: Main menu "Calibration"

▷ At initial start-up of the system, if the function is enabled for the first time, the functional elements can show the state "Off" or "Measure". After calibration they will be switched over to either "Totalize" or "Measure".

5.2 Operation of the Calibration Function in General

- ▷ The main function "Calibration" is enabled with the main function key in the bottom section. The submenus, e.g. for the calibration of sensors or totalizers resp. supplies of media (pumps) are accessible from the main display.

You can enable the function and change calibration parameters, such as for the re-calibration of the pH-sensor or changing the totalizer parameters, at any time as required even during a process.

- ▷ All inputs for a calibration function are done in a logical order of steps via submenus, which are either displayed automatically or need to be enabled via the mode selection key "Mode".

The TAG of the function key "Mode" switches from "meas" to "calib" after starting the first step of the calibration routine and returns to "meas" after the last step. The color changes from grey (state "off" or "inactive") to red ("Calibrate") and finally to green ("Measure"):

Fig. 5-2: Key display during calibration of sensors

- ▷ In some calibration menus context based information guides through the routines, telling which inputs and selections are to be made. The system recalculates the calibration-dependent (measured or calculated) process values after completion of the routine, then.
- ▷ All calibration parameters will be stored permanently. When the control unit or controlled device are switched-off and restarted for continuing a process or for a new process afterwards, the system uses these parameters again for the calculation of process values, until the calibration is done again.

5.3 pH Calibration (Conventional pH sensors)

The conventional pH sensors are calibrated in a 2 point calibration. This determines the sensor parameters "zero deviation" and "slope".

The calibration function implemented in the DCU system for conventional pH sensors allows that you can

- ▷ enter the zero pH and voltage directly into the system
- ▷ enter the slope pH and voltage directly into the system
- ▷ proceed through a complete calibration course
- ▷ recalibrate the pH-sensor at any time during the process
- ▷ do either a zero calibration or a slope calibration separately.

During pH measurement, the system calculates the pH value according to the [► Nernst equation] from the sensor potential, considering the zero deviation, slope and temperature. During the calibration you can enter the reference temperature manually, during the pH measurement the temperature compensation is carried out automatically, based on the temperature measured in the medium.

To re-calibrate the pH sensor, you can measure the pH value in a sample taken from the process and enter the value into the calibration menu. The calibration function compares the pH value measured online with the one determined externally, calculates the resulting zero offset and displays the corrected process value.

In addition to the pH value, the sensor's operator screen also displays the sensor voltage and the sensor parameters "zero offset" and "slope". This allows you to easily check the integrity of the pH sensor.

5.3.1 Overview on the pH Calibration Function


Fig. 5-1: Submenus of the pH calibration function

Tag	Entry	Function, display, required input
Measure		Switch to pH measurement after calibration routine is completed
Calibrate		Starting the calibration routine
Re-Calibrate		Starting the re-calibration
Calibrate Zero		Zero point calibration as a single step
Calibrate Slope		Slope calibration as a single step
pH	pH	pH value and or entering the pH of the external sample during re-calibration
Sensor	mV	Combination sensor voltage (raw signal)
TEMP	°C	Setting for temperature compensation
Zero	mV	Display the zero offset
Slope	mV/pH	Display the slope
Manual		Manual mode
Auto		Auto mode

5.3.2 Proceeding through a Calibration

For proceeding through a complete calibration routine including the "zero" and "slope" calibration please note the steps a), b) and c) of fig. "Submenus of the pH calibration function".

- ▷ To run the calibration "manual" temperature compensation select "Manual" at step a) and enter the temperature which the system should apply for the calibration routine.

With "Auto"matic temperature compensation the system uses the temperature measured at the test point assigned.

- ▷ If only the sensor's "Zero" should be calibrated, select "Calibrate Zero" and observe step "b)". This may be necessary if there is a longer break between "Zero" and "Slope" calibration, while other operations need be performed at the controlled device.
- ▷ To calibration the "Slope" only, select "Calibrate Slope" and observe step "c)".

5.3.3 Re-calibration

The re-calibration functions allows to adapt the calibration to a changed measuring behavior of the pH sensors effected by impacts of the medium during the process, as needed.

1. Measure the pH value of a sample taken from medium (container). Use a precise and carefully calibrated external measurement device.
2. Tap **[Measure]** and select **[Re-Calibrate]**. Enter the pH measured in the sample:
 - ▷ The DCU system calculates the resulting zero offset and displays the corrected pH value.

5.3.4 Special Notes

- ▷ If possible, use buffer solutions available from the manufacturer of the sensor, as supplied with the pH sensor.
- ▷ If "zero offset" and "slope" values are known, you can enter them directly into the relevant fields.
- The lifetime of the pH sensor depends on the process conditions and is limited. A pH sensor should be serviced and replaced if these values exceed the specified range¹:
 - zero offset beyond the range of -30 ... +30 mV or
 - slope beyond the range of 56 ... 59 mV/pH.

¹ Note the [▶ manufacturer's information] for valid limit thresholds.

6. Main function "Control Loops"

6.1 Overview on Controllers

6. Main function "Control Loops"

The controllers included in a system depend on the controlled device and the configuration ordered by the customer.

A system configuration considers controllers and respective settings suitable for the controlled device [► configuration documentation].
A controller has such control parameters preconfigured which ensure a stable control behavior.


If required by the process, single control parameters or the control structure (e.g. PID, if available) can be adapted to the control task.

Since the FlexAct® System has only the following control function there is no extra "Main Function" implemented; the control settings can be run from the "Main" screen.

Controller	Type of controller, actuator, control function
Pump controller	Setpoint controller for triggering an external dosing pump

6.1.1 Access to Controller Menus

The FlexAct® pump control function is accessed directly from the "Main" menu, e.g. for input of setpoints or the selection of the controller modes.

1. After tapping of  in the menu "Main" a submenu for input of the setpoint and the selection of the operation mode appears.
2. Tapping of  opens the submenu for setting the specific control parameters, such as the alarm monitoring, and the parameterization of the pump controller [► section "Pump controller"].

► Example: Access from „Main" menu to the pump controller:

Fig. 6-1: Pump control menu accessed from main menu

Controller functions which are accessible for routine operation by normal users include the selection of the controller modes as well as the input and change of setpoints and of controller outputs.

Inputs of "Setpoints" must be within the limits of the controlled variable. Invalid entries (which exceed the limits) will be ignored.

During "Remote" operation the host computer transmits the setpoints and operation modes.

Settings not required for routine operation are entered via parameterization functions. If available in the configuration this includes the setting of the controller output limits "Min"/"Max", or, in case of PID controllers, of a deadband DEADB and of the PID parameters [► section "Controller Parameterization in General"].

The controllers switch over the modes of their controller actuators at bumpless transfer.

Tags and entries of the menus (Overview)

Not all tags and entries are available for all controllers.

Tag	Entry	Function, display, required input
Mode		Controller modes
	off	► the controller (and assigned servo controllers) are switched off; the output is in its defined state (stop state of actuator)
	auto	► controller is enabled and triggers the actuators according to the preset setpoint (external setpoint at "Remote" operation) – setpoints are entered via the controller menu entry "Setpoint"; – servo controllers are switched over to "cascade" mode
	manual	► manual access to controller output (manual setting of output)
	case	► controller is servo controller in control cascade, the master controller delivers the setpoint (signal) – no setpoint can be entered directly for this servo controller – mode "cascade" is not possible for all controllers
	profile	► controller operates a timer profile – mode "profile" is not possible for all controllers
Prozess value	[PV]	Current process value in its assigned unit, e.g. physical unit
Setpoint	[PV]	Setpoint of process value in its assigned unit, e.g. physical unit
OUT	%	Controller output; inputs are only possible in mode "manual"
Cascade	[controller]	Servo controller triggered by the output of the master controller
		Access to parameterization menu of the controller, usually via the standard password [► supplement]; not visible while entered
ALARM		Enabling of the alarm display
PROFILE		Input menu for a setpoint profile (if available in the configuration)

Unchecked, unauthorized changes of control parameters can cause instable controls and unexpected erroneous function of actuators.

Changing of controller parameters is restricted to authorized persons only. The parameterization displays are accessible only by entering a password. The password must only be forwarded to persons which are aware of the effects of changes of the control parameters.

Usually controller parameters are not changed. An exception are parameters effected by the process, such as of pH and pO₂-control.

In all configurations the controllers are parameterized so that they (under normal conditions) can provide a stable behavior of the controls of the devices they are assigned to. Parameters preset in the factory are shown in the [► configuration documentation] of the delivered device.

To adapt a controller to the behavior of its control loop it may be possible to change these characteristic values in the parameterization display:

- ▷ MIN, MAX Min. and max. controller output limits
- ▷ DEADB Setting for deadband (only available for PID-controllers)
- ▷ XP, TI, TD PID parameters (only available for PID-controllers)
- ▷ Parameterization menus are accessible via the touch-key " " in the controller menu and input of a password [► Supplement].

6.3.1 Tags and Entries of Parameterization Displays (Overview)

Not all tags and functions are available for all controllers.

Tag	Entry	Function, display, required input
MIN (Minimum) MAX (Maximum)	 as "%"	Minimum and maximum output limits; in cascade control: min/max of servo controller, switching over condition for previous next controller.
DEADB	 unit "PV"	Deadband; the control remains inactive as long as the current process value deviates from the setpoint for less then DEADB.
Xp	 as "%"	P-portion (proportional range); the amplification of the output signal is proportional to the input signal (deviation of current value from setpoint).
TI	 as "%"	Integral portion; time aspect of the controller action; the controller action slows down at increasing I-portion (and vice versa).
TD	w.v.	Differential portion; "damping" of the control action; the controller action is weakened at increasing D-portion (and vice versa).
TD	s	Rate (delay) time
		Controller assigned in cascade, order according to priority – [Tag] = short name of the controller or actuator
Hyst. time	mm:ss	Switching delay time (characteristic value) of the controller
Process Value		Selected process value which should be controlled (if more than one)

6.3.2 Output Limits

Safe operation of the controlled device requires that the output of some controllers must be limited [► Configuration documentation].

Limiting the controller outputs avoids unintended actuator changes and, in case of cascade controllers, limits the setpoint for the servo controller.

- Limit settings are made in terms of “%” of the max. controller out.
- For full controller output range the following limits must be set:
 - single controller outputs: MIN = 0%, MAX = 100%
 - splitrage controller outputs: MIN = -100%, MAX = 100%

6.3.3 Deadband

For current values which vary stochastically a deadband provides a stable control at minimum corrector actions. In case of splitrage controller outputs this avoids an oscillation of the outputs. As long as the control deviation is within the deadband the controller output remains constant.

- The deadband is set in the entry DEADB in the dimensions of the process value or setpoint. It is symmetrical to the adjusted setpoint.
 - Example for pH controller: adjusted setpoint 6.0 pH, measurement range pH = 4 ... 10 pH; measurement span 6 pH
 - adjusted deadband 0.3 pH.The pH-control remains inactive within pH 5.7 - 6.3.

6.3.4 Control Structure and Optimization of PID-Controllers

The control structure results from setting certain PID-parameters to zero:

P controller	TI = 0, TD = 0
PI controller	TD = 0
PD controller	TI = 0
PID controller	all PID parameters specified

The setting of PID-controllers can be optimized via their PID-parameters [► above table] The optimization requires knowledge in control theory [► technical literature, e.g. rules by Ziegler Nichols]. Guidelines are:

- Enable the D section (TD) only at relatively stable current values. If current values vary stochastically TD changes the controller output considerably and fast. This will cause an unstable control behavior.
- As a general rule the “TI:TD” ratio should be 4:1.
- In the case of a periodic oscillating control loop you can increase both Xp and the ratio “TI:TD”.
- If a current value reaches the setpoint too slow after setpoint jumps or if it drifts despite of the control, you can decrease Xp or “TI:TD”.

6.4 Pump Controller

The pump controller can trigger a continuous or discontinuous pump for transfer of the medium of substrate via a corresponding setpoint signal.

- ▷ At the FlexAct® systems the pump controller triggers an external pump which feeds the medium from the source (e.g. LevMixer® Pallettank®) via the filters into the receiving Pallettanks® for storage. The volumes to be transferred can be specified via the weight "FLT-WGHT" of the receiving Pallettanks®. The amount will be calculated from the weight reduction of the source tank.
- ▷ The filtrate transfer is provided as phase control. It can be repeated as often as possible, until the source tank has a minimum weight.

6.4.1 Submenus and Operating Displays




Fig. 6-2: Overview on the pump control function

6.4.2 Settings for Pump Operation

Tag	Entry	Function, display, required input
Operating display		
Mode		Input (selection) of controller mode
		– controller switched off, pump off
		– controller switched on, automatic operation towards setpoint
		– manual access to controller out
Pump	%	Current (signal of) pump delivery, as percentage value of max. delivery
Setpoint	%	Controller setpoint for the pump
Out	%	Current controller out signal as percentage value of max. out
		Access to parameterization menu via standard password [► supplement]
Alarm Pump		Monitoring of the course of the current values vs. the alarm limits
Hi Lim	%	– high limit alarm
Lo Lim	%	– low limit alarm
Alarm	state	– state of alarm monitoring, switching over "disabled" – "enabled"
Parameterization Display		
Lower limit	%	Lower speed limit, which shall not be reduced any more
Upper limit	%	Upper speed limit, which shall not be exceeded.

Notes on pump operation:

- ▷ The pump can only be switched to and run in manual mode, if no phase [► Filtration] is active.
- ▷ The pump can be switched to run in reverse direction via the  key in the main menu.

6.4.3 Parameterization of the Pump Control

The parameterization display can be used to limit the output to a certain minimum or maximum pump speed.

- ▷ The MIN | MAX-output limits are entered as percentage of the control range of the pump speed.

6.4.4 Special Notes

The pump will be stopped by a corresponding interlock function if:

- ▷ the contents of the source container falls below its preset minimum
 - the source container is monitored by an external scale which transmits the corresponding weight signal to the DCU-system
 - the alarm limits are set in the DCU alarm display for the scale.
- ▷ if the filling volume of the filtrate bag, exceeds the set value
 - the filling volume is set via the "SET_WGHT" parameter in the [► Parameterization display of the filtrate phase].

7. Phases, Step Controls

7.1 Overview on Phases

"Phases" provide access to combined functions for device or process controls by predefined sequences of process steps. These sequences can run automatically or be combined with manual operations, which the operator has to perform and confirm at the DCU system.

7.2 Phases of the FlexAct® System

In the current configuration of the DCU for FlexAct® systems phases are not implemented as extra "Main function", they are accessible from the Main menu, instead. The configuration provides a single phase for buffer filtration, consisting of two subphases:

- ▷ Initial filling of the system [▷ "Tube filling"]
- ▷ Filtration of a certain buffer volume [▷ "buffer filtration"], which can be repeated as often as possible from the LevMixer® Pallettank®.


Fig. 7-3: Access to submenu "Phase Filtration" from "Main" menu

The DCU executes the phase semiautomatically, i.e. the operator will start the phase and then confirm the fulfillment of predefined steps, in order to branch to an alternate subphase or to proceed with the next step.

The operator can start and stop a phase at any time. He can run a phase step-by-step. Each step represents a control action which triggers the process equipment assigned, e.g. changing a setpoint of a controller.

7.2.1 Overview on the Filtration Phase


Fig. 7-4: Displays with steps and transitions for running the phase and parameterization

Tags and entries of the menus

The displays, tags and entries shown here are only exemplary. For a complete list refer to the [► configuration documentation, section "Sequences"] supplied with the system.

Tag	Entry	Function, display, required input
Operation displays (running the phase):		
State		Current mode of the phase:
	idle	– idle: phase inactive
	running	– running: processing the phase (all steps)
Runtime	hh:mm:ss	duration of step after start
actual step		Explanation of current step; request for confirmation of step, e.g.:
condition	Q_install	– verifying, whether installation is complete and correct;
next step	Q_first use	– decision of whether to branch to tube filling or filtration steps;
	Q_filling	– verifying, whether the setting for filling is correct, starting the step;
	S_filling	– system will be filled, stop by the operator after completion;
	Q_filtration	– verifying, whether the setting for filtration is correct;
	S_filtration	– buffer will be transferred (filtered);
	S_slow filtr	– slow down of filtration to allow for stop in order to prevent excess filling of bags, prior to release an overfilling cut-off;
	Q_repeat	– request on running another filtration or ending the phase;
	S_End filtr	– branch to ending or repeating the filtration phase.
	<input type="checkbox"/> Yes / <input type="checkbox"/> No	– Branching to alternate subphases as shown in "Condition" section.
FLT_WHGT	kg	Current weight of filtrate bag calculated from weight reduction of the source (e.g. LevMixer® Palletank®).
Parameterization displays (accessible only by entering the "Standard password" [► Supplement]):		
Set_WGHT	kg	Set weight for filtrate bag.
SET_Pump	%	Set value for pump delivery (percentage of delivery range).
S_Pump	%	Set value for pump at the "S-filling" step; provides that the filling can be stopped precisely after filtration of the intended volume; prevents excess filling of the filtrate bags.

During the filtration phase related “Condition” notes inform about checking the parameters of the pump, the LevMixer® Palletank® and the receiving bags and request a verification whether all equipment is properly connected [► Operating Manual “FlexAct® Central Operating Module COM”]. The preparatory measures are described in details in the [► Process description “Buffer preparation”]. In general, however, it is recommended to do all checks in advance of running the phase.

Prior to running the filtration phase you should do the following:

- ▷ Check whether “SET_WGHT” on the “Main” display corresponds with the filling volume intended for the filtrate bags.
- ▷ If not, enable the parameterization display and set the “Filtration SET_WGHT” as required.
- ▷ Check in the parameterization display whether Filtration SET_PUMP and Filtration S_PUMP are set as intended or change, as necessary.

To prevent from unauthorized changes, the parameterization display is password-protected [► section “Password” in the supplement].

The filtration phase entails steps for preparation of the assembly, which includes the prefilling of the filters, and the steps for the filtration of the buffer from the LevMixer® Palletank® into the receiving bags. During the step “Q_First Use” the operator can branch to the steps for filling the tubes or the steps of buffer filtration.

After completion of the filling steps the phase will automatically switch over to the first filtration step and wait for a confirmation to start.

The system allows a precise filling and manual stopping of the transfer, if the intended volume is introduced into, e.g., the receiving bags. This is done by reducing the pump speed from the “SET_PUMP” setting to the “S_PUMP” setting automatically when the amount of buffer removed from the LevMixer® Palletank® reaches the setpoint “SET_WGHT”.

After running the filtration once, the operator can branch to repeating the filtration with another filtrate bag or to end the filtration at all.

The system monitors the filtration against overfilling of the filtrate bags and emptying of the LevMixer® palletank below a preset minimum volume. Corresponding interlocks will cut-off the pump if

- ▷ the weight reduction of the LevMixer® Palletank® exceeds the weight which is set in the parameterization display of the filtration phase as SET_WGHT, or
- ▷ if the weight of the LevMixer® Palletank® falls below the “WEIGHT Low Lim” alarm limit set in the corresponding “Alarm” display [► Setting and enabling the weight alarm monitoring].

8. Main Function "Settings"

This main function permits changes of the system configuration. Changes of settings, which are unsuitable or not permissible for a controlled device, can cause unpredictable malfunctions.

Settings which effect the safe operation of the controlled device, are password-protected to prevent unauthorized access. The standard password [► Supplement] must only be passed on to authorized users. The service password [► supplied extra] must only be passed on to service personnel or authorized users.

8.1 Overview

The DCU system provides several functions for system maintenance and trouble shooting with the main function "Settings":

- General settings like date, time, the recovery time "Failtime" after errors, enabling of a screen saver or parameterization of the communication with external devices ["Internet Configuration"].
- Setting of process values "PV" and their ranges and limits.
- Manual operation of in- and outputs for simulation.
- Service functions, such as a system reset or for selection of a different system configuration, as far as multiple configurations are implemented and accessible for use.


Fig. 8-1: Main menu "Settings"

Tag	Entry	Displayed value or function, required input
Hardware	Microbox XXXXY	Hardware of the DCU system, e.g. Microbox PC420
Firmware	X.YY	Firmware version implemented in the system
Configuration	Version_ Unit_ Config	Configuration implemented in the system

If you require further information or in case of a claim on malfunctions please refer to the firmware and configuration of your system.

8.1.1 Settings Display and Submenus of DCU System for FlexAct®


Fig. 8-2: Display and submenus of main function "Settings"

- Some settings will be set once at initial start-up or after a system reset; the related submenus will not be accessible later on.
- Some settings (the access to related submenus) are restricted to authorized users only.

8.2 Setting of
System Parameters

8.2.1 General System
Settings

The "Settings" function "System Parameters" can be used for general settings of the DCU system.

8.2.2 Operating Displays

The "Time" and "Date" can only be changed within 5 minutes after start-up of the system.

For details see the sequence of displays shown in the fig. "Display and submenus" of the main function "Settings" in section 8.1.

Tag	Entry	Displayed value or function, required input
Time	hh:mm:ss	Display of current time, format is [hours:minutes:seconds]
Date	dd.mm.yyyy	Display of current date, format is [day.month.year]
Failtime	hh:mm:ss	Input of power failure time (Failtime) for the system behavior after a restart. The system will behave as follows: <ul style="list-style-type: none">– power failure < Failtime: system continuous with the last settings– power failure > Failtime: systems returns to default settings
Screen saver	hh:mm:ss	Input of activation time of a screen saver. Access to the display will only be possible by entering the standard password.
Internet configuration		Configuration of the Ethernet interface for host connection by means of the LAN applied in the laboratory, input of
	IP-Address	– 12-digit IP-address, according to the local LAN specifications
	Gateway	– 12-digit Gateway-address, according to the local LAN specifications
	Netmask	– Netmask 12-digit Netmask-address, according to the local LAN specifications

8.3 Setting of Process Value Ranges

Only authorized personnel shall change the settings in this menu.
The settings effect the setpoint limits of the controllers, for instance.


Fig. 8-3: Overview on setting of process value ranges (examples for selecting submenus)

- ▷ The PV ranges configured for the delivery state of the controlled device are specified in the [► "Configuration documentation"]. The "Min" and "Max" settings of all process values, the alarm limits and the alarm monitoring mode can be changed as required within the limits of the configuration.

8.4 Manual Operation

For testing and trouble shooting of the DCU system and of the controlled device the analog and digital inputs and outputs, the controllers and digital counters can be switched over to manual operation. This disconnects the DCU functions (e.g the inputs of measured signals, outputs of control signals) from their devices and allows for direct input of parameters and control signals as suitable for simulation of the function.

Settings entered in manual mode have top priority, they will override all other settings and functions of inputs and outputs of the DCU system.

8.4.1 Selection of the Operating Displays (Overview)


Fig. 8-4: Overview of the displays "Manual Operation"

8.4.2 Manual Operation of Digital Inputs

For manual operation disconnect the digital input from the external device, such as a limiter. You can simulate the input signal via entering "on" or "off" for the state of the input.

Operating Displays



Fig. 8-5: Manual operation of digital inputs
Switching to manual mode and returning to normal operation (example)

Menu Tags and Entries:

Tag	Entry	Function, required input
IO Tag	FILL_STOP, 1DInn, etc	Display of TAG of digital input from the process – FILL_STOP: Input assigned, e.g. the status signal from a second remote scale, in order to allow for an High Level Alarm. – 1DInn: no input signal assigned to channel no. "nn".
Mode	auto manual	Current mode of digital input channel: – normal operation, external input is transmitted to DCU – manual operation, digital input is set manually for simulation
Value	off on	Specified input and state (green : automatic; yellow : manually operated) – digital input switched off – digital input switched on

Operating Information

1. If the table has more entries than visible in one display scroll down with the slider to the entry which should be changed. Tap the cell in the column "Value", to enable the submenu for selection of mode.
2. Tap "Manual" and enable the touch key of the intended mode of the digital input.
3. Confirm with . The entry below "Value" changes correspondingly. The color of the cell changes to yellow to indicate, that the input is operated manually for simulation.

TAGs for digital process inputs

- ▷ A complete list of all TAGs is included in the supplement. Unassigned reserve channels are defined by the generic name of the channel, with the list no. of the channel.

Special Notes

- ▷ For the state of the inputs the following signal levels are applied:
- OFF : 0 V
 - ON : 24 V für process inputs (DIP)

After operation in "manual" mode all related digital inputs must be reset to "auto". Otherwise the normal functions of the DCU system will be restricted. For proceeding see above figure.

8.4.3 Manual Operation
for Digital Outputs

For manual operation you can disconnect a digital output from its internal function and adjust manually. You can switch static outputs, such as the heater controller, "on" and "off" [► fig. below, course on right hand]. For pulse-width modulated outputs you will enter the pulse/pause relation per operation cycle in terms of "%" [► course on left hand].

Operating Displays


Fig. 8-6: Manual operation of digital outputs;
Switching to manual mode and returning to normal operation (example).

Menu Tags and Entries:

Tag	Entry	Meaning
I O Addr.	PUMP_AUTO, 1DOnn, etc.	Display of TAG of digital output to the process – PUMP_AUTO: Output assigned, example: status signal to pump – DOPnn: no output signal assigned to channel no. "nn"
Mode	auto lock shut manual host opst	Current mode of digital output channel: – normal operation, DCU triggers the output – (display only) output is locked, such as due to alarm locking – (display only) output in safe state via shut-down function – manual operation, manual setting of digital output – output triggered by host computer – output in default state [► configuration documentation]
Value	off on nn%	Specified input and state (green: automatic; yellow: manually operated) – digital input switched off – digital input switched on – pulse-pause relation for pwm output, in the range of 0 ... 100%

Operating Information

1. Scroll down with the slider to the entry which should be changed. Tap the table entry below "Value" to open the submenu.
2. Select the "Manual" mode. For static outputs enable the intended mode. For "pwm outputs" enable the touch key of the process value. The submenu for data input opens.
3. Confirm the selection or input with **[ok]** changes the "Val" entry accordingly. The color of the table cell changes to yellow.

Special Notes

- Several functions can trigger one digital output. The active function is displayed in the "Mode:..."-section in this order of priority:

Highest priority	Lowest priority
– shut Down	– operating state (opst)
– manual operation	
– locking	
– pump calibration	
– controller, timers, sensors, scales	

- The status can have the following signal levels:
- off: 0 V
 - on: 24 V for process outputs (DOP)
- For pulse-width modulated (pwm) outputs the pulse control factor is displayed or set. The cycle time is set in the configuration:

At a cycle time of 10 sec. and a pwm output of 40%, the output will be "on" for 4 s and "off" for 6 s.

After operation in "manual" mode all related digital outputs must be reset to "auto". Otherwise the normal functions of the DCU system will be restricted. For proceeding see above figure.

8.4.4 Manual Operation of Analog Inputs

For manual operation disconnect the analog input from the external device, e.g. an amplifier. You can simulate the input via entering a voltage signal of 0 ... 10 V or a current signal of 4 ... 20 mA, corresponding to a relative signal level of 0 ... 100% of the physical range of the input.

Operating Display


Fig. 8-7: Manual Operation of Analog Inputs
Switching to manual mode and returning to normal operation (example).

Menu Tags and Entries:

Tag	Entry	Function, required input
IO Tag	pH, 1Alnn, etc.	Display of TAG of analog input from process. Examples are: <ul style="list-style-type: none"> – pH: input assigned, signal from pH measurement – Alnn: no input signal (channel) assigned to analog input "nn"
Mode	auto manual	Current mode of analog input <ul style="list-style-type: none"> – normal operation, input from device triggers DCU – manual operation, setting of analog input for simulation
Value	0 ... 10 V 4 ... 20 mA	Specified input and state (green: automatic; yellow: manually operated) <ul style="list-style-type: none"> – physical signal representing a relative input level within 0 ... 100%

TAGs of Analog Inputs

- An overview of tags of assigned inputs is provided in the supplement. Unassigned reserve channels are defined by the generic name of the channel, with the list no. of the channel.

Operating Information

1. If the table has more entries than visible in one display scroll down with the slider to the entry which should be changed.
2. Tap the table entry below "Val" to open the submenu for selection of the mode.
3. Select the "Manual" mode. The submenu for data input opens at tapping the parameter key.
4. Enter the value for the analog input and confirm with **[ok]**.
The "Val" entry is changed accordingly, the table cell is displayed yellow, indicating that the input is controlled manually.

Special Notes

- Analog inputs (AIP) can be simulated by voltage or current signals within the following signal levels:
 - 0 ... 10 V (0 ... 100%)
 - 0 ... 20 mA (0 ... 100%)
 - 4 ... 20 mA (0 ... 100%)
- To which physical process values this corresponds depends on the range of its measured signal. For simulation of a specific physical input signal of an analog input, the related voltage or current signal level within 0...100% must be evaluated.

After operation in the "manual" mode all related analog inputs must be reset to "auto". Otherwise the functions of the DCU system will be restricted. For proceeding see above figure.

8.4.5 Manual Operation
for Analog Outputs

You can disconnect all analog outputs from the internal DCU functions and adjust the output signal manually within a physical voltage output signal in the range of 0...10 V or a current output in the range of 4...20 mA, which corresponds with a relative signal level of 0 ... 100%.

Operating Display


Fig. 8-8: Manual Operation of analog outputs
Switching to manual mode and returning to normal operation (example).

Menu Tags and Entries:

Tag	Entry	Function, required input
I O, Tag	PUMP, 1AOnn, etc.	Tag of analog output, manual input of mode and output signal – PUMP: output assigned, output signal for extern pump – AOPnn: no output assigned to analog output channel "nn"
Mode	auto lock (display only) shut (display only) manual	Display of current mode – normal operation, DCU triggers analog output (remote unit) – analog output switched to safe state via locking function – analog output switched to safe state via shut down function – manual operation with signal as specified for simulation
Value	0 ... 10 V 4 ... 20 mA	Specified input and state (green: automatic; yellow: manually operated) – physical signal representing an output of 0 ... 100%

Operating Information

1. If the table has more entries than visible in one display scroll down with the slider to the entry which should be changed.
2. Tap the table entry below "Val" to open the submenu for selection of the mode.
3. Select the "Manual" mode. The submenu for data input opens at tapping the parameter key.
4. Enter the value for the digital output and confirm with [ok]. The "Val" entry is changed accordingly, the table cell is displayed yellow, indicating that the output is controlled manually.

Special Notes

Outputs have the following order of priority:

- | | |
|---|---|
| <p>▷ Highest priority:</p> <ul style="list-style-type: none"> – shut down – manual operation – locking | <p>▷ Lowest priority:</p> <ul style="list-style-type: none"> – controller action (setting), etc. |
|---|---|
-
- ▷ A list of tags of analog outputs is provided in the supplement. Unassigned reserve channels are defined by the generic name of the channel, with the list no. of the channel.
 - ▷ Physical signal levels of analog outputs can be configured within:
 - 0 ... 10 V (0 ... 100%)
 - 4 ... 20 mA (0 ... 100%)
 - ▷ To which physical process values (e.g. stirrer speed) this corresponds depends on the range of its measured signal. For simulation of a specific physical output signal of an analog output, the related voltage or current signal level within 0...100% must be evaluated.

After operation in "manual" mode all related digital outputs must be reset to "auto". Otherwise the functions of the DCU system will be restricted. For proceeding see the above figure.

8.4.6 Manual Operation of Controls

For simulation the operator can disconnect the controllers from their DCU control functions. He can simulate the controller signal via entering a variable within the "Min" and "Max" limits of the process value, as permitted and provided by the configuration.

Operating Displays:


Fig. 8-9: Manual operation of controls for simulation

Menu Tags and Entries:

Tag	Entry	Function, required input
CV	PUMP, etc.	Display of TAG of control function controlled variable
Mode	[off] [auto] [manual]	Current mode of the control function – control function disabled – normal operation, DCU triggers control function – manual operation, input of controlled value for simulation
Value	[PV]	Signal specified for simulation – variable (process value, setpoint or controller output), in its physical dimension

Operating Information

1. If the table has more entries than visible in one display scroll down with the slider to the entry which should be changed.
2. Tap the table entry below "Val" to open the submenu for entering the PV for simulation.
3. Select the "Manual" mode. The submenu for data input opens at tapping the parameter key.
4. Enter the value for the digital output and confirm with [ok]. The "Val" entry is changed accordingly, the table cell is displayed yellow, indicating that the control has a manual setting for simulation.

Special Notes

- ▷ For controlled variables the signal level can only be configured within the range of the "Min" and "Max" limits of the process value, see the displays and the [▶ configuration documentation].

After operation in the "manual" mode all related analog inputs must be reset to "auto". Otherwise the normal functions of the DCU system will be restricted. For proceeding see the above figure.

The service menu is only accessible via a special service password. The password must not be made accessible to unauthorized personnel.

Rebuilding of internal modules of the control unit and the update of the configuration for operation beyond the original specifications must only be done by authorized service personnel or in the factory.

A "System Reset" and loading of a new configuration deletes all user-modified settings and enables the default delivery state. Therefore you should do a system reset only if absolutely necessary, such as for clearing of malfunctions, which cannot be solved otherwise.

After a system reset verify whether the specific settings required for the controlled device, are properly set, such as the output limits of controllers [► "Configuration documentation"].

- The "Service" menu provides access to the configurations delivered with the DCU system and allows for a system reset.
- Implemented configurations include parameter sets which ensure the proper operation of the controlled devices they are assigned to [► Configuration Documentation]. For customized units the extent of functions can be specified in the [► Functional Specification].
- A changed or new configuration may be required in these cases:
 - at operation beyond the original specifications;
 - at use of the DCU-system with different controlled device(s).

System-Reset and Changes of a Configuration

1. Prior to a system reset or loading of another configuration please make a record of the user-defined settings, which may be required for reproducing or re-entering the parameters of your process, e.g.:
 - measurement ranges,
 - controller parameters,
 - process value and digital input alarms,
 - communication parameters.
2. Enable the touch key "Service" in the "Settings" menu. Enter the service password. The current configuration is displayed.
3. Tap the touch key "Configuration". A list of all configurations available in the system appears. Select the configuration.
4. If the intended configuration is not displayed in the list, please contact the service representative of Sartorius Stedim Systems GmbH.
5. Confirm the loading of the new configuration with [ok].
 - reloading of a current configuration has the same effect as a "System Reset". It therefore reloads only the default settings. User-defined settings need be reentered from the records.
 - if no new configuration should be applied, you can quit the menu without changes either via key "C" or the button [X].
6. Re-enter the user-defined settings of the old configuration, as far as required for the intended process, into the corresponding menus.

9. Supplement

9.1 Alarms and Messages

9. Supplement

If an alarm occurs, a corresponding message is displayed on-top of the current display and in the headline. The alarm text in the headline and the alarm symbol in the main function key "Alarm" change their background color to red [🔔 → 🔴🔔]:


Fig. 9-1: Alarm message at occurrence of an alarm "Factory reset"

The main function key "[Alarm 🔴🔔]" provides access to a table showing all alarms which occurred during the process:



Fig. 9-2: List of all alarms which occurred, with a short description

9.1.1 Operating Information

- The operator can handle alarm messages as follows:
- 1. Confirmation of the alarm in the pop-up window with "Acknowledge". The alarm display is closed and the red background of both the alarm text and symbol disappears. The alarm is listed as "ACK" (acknowledged alarm) in the table of alarms.
 - 2. Closing the pop-up window via  without acknowledgement. Alarm message and symbol are kept highlighted with the red background. The alarm is listed as "UNACK" in the table of alarms.
 - 3. Opening the alarms table via the main function key "". Unacknowledged alarms can be acknowledged with "ACK". Tapping "RST" ("Reset") deletes the message from the list of alarms.
- ▷ How to handle alarms and other messages depend on their cause and meaning. The tables below give an overview on this. If questions arise on alarms and appropriate remedial measures, please contact your representative of Sartorius Stedim Systems GmbH.

9.1.2 Overview on Process Alarms

Available alarms and assigned alarm messages depend on the controlled device and on the specific configuration. For detailed information refer to the [\[► Configuration documentation\]](#)

Text of alarm line	Meaning	Remedial measures
"Name State alarm ..."	alarm of digital input "name"	confirm alarm with "ACK"
"Name low alarm ..."	the process value "name" is below its lower limit threshold	confirm alarm with "ACK"
"Name high alarm ..."	the process value "name" is below its upper limit threshold	confirm alarm with "ACK"
"Fuse blown ..."	fuse blown or automatic cut-out is enabled	replace fuse, switch on automatic cut-out again; if alarm is repeated contact customers service
"TEMP Sens. failure ..."	Pt-100 or cable is defective	check cable plug connection; if in good order: replace sensor

9.1.3 Process Messages

Available process messages depend on the controlled device and on the specific configuration. For detailed information refer to the [► Configuration documentation]

Text	Meaning	Measures
"Shut down ... "name" ..."	Emergency cut off – key of controlled device enabled	Switch-on again the controlled device "name" via the "shut down" key or the Smains switch
"Name <3mA alarm ..." (name: TAG of analog input)	The analog input "name" has a broken cable	Check cable connection and peripheral unit connected
PUMP_FAIL	Pump not ready, e.g. pump head lid open, tubing not inserted properly	Check pump, whether the tubing is inserted properly and the pump head lid is closed.

9.2 Trouble Shooting for the DCU system

If problems arise with the use of the DCU system (malfunctions, technical problems), Sartorius Stedim Systems GmbH needs a detailed description of the problem.

The operator or the customers' service can use the form "Function bug handling" which is available on request, to describe the problem. This will enable Sartorius Stedim Systems GmbH to reply within the shortest possible time offering a solution for the problem.

In case of problems we recommend that the user or his service should contact the responsible local service representative of Sartorius Stedim Systems GmbH. However, if no local service representative is available, he can directly send the completed form to:

- Sartorius Stedim Systems GmbH
Dept. FA (Product engineering, automation projects)
Robert-Bosch-Straße 5 – 7

D-34302 Guxhagen, Fed. Republic of Germany
phone +49.5665.407.0, fax +49.5665.407.2200

9.3 Locking Functions

- Locking functions cannot be changed by users. For any output which is effected by a locking function the system displays the state "lock" at manual operation.
- The extent of locks is specific for the controlled device and defined in the configuration. For details refer to the [► Configuration documentation], which can either be supplied with the control unit, be included in the [► Technical Documentation] or be supplied on request, respectively.

9.4 Abbreviations

A DCU system does not necessarily include all TAGs and entries shown here. This depends on the configuration [► Configuration documentation].

9.4.1 Process Values

Tag	Meaning
TEMP	temperature, signal from temperature sensor, e.g. Pt-100
pH	pH value, signal from pH-electrode
PRESS	pressure, signal from pressure sensor
COND	conductivity, signal from conductivity sensor

9.4.2 Digital Inputs

Tag	Displayed value or function, required input
PUMP_FAIL	failure message from pump; requires check of pump and cable
FILL_STOP	filling of filtrate bag stopped by external (scale) signal, e.g. if weight of filtrate bag exceeds the setpoint or if weight of supply container (media source) falls below its "Low alarm" limit.

9.4.3 Digital Outputs

Tag	Displayed value or function, required input
PUMP	signal for external pump, e.g.: <ul style="list-style-type: none">– PUMP_AUTO: Pump running as set on the pump– PUMP_STOP: Pump interlock – stop if weight of filtrate bag exceeds the setpoint or the media container is emptied below its minimum weight (min. alarm limit)– PUMP_REVERSE: Pump running in reverse mode

9.4.4 Analog Outputs

Tag	Displayed value or function, required input
PUMP	Controlled output signal for external pump, as percentage of the pump speed range

9.5 Arrangement and Assignment Plans

- ▷ For the arrangement and assignment plans of the control unit please refer to the "Technical Documentation DCU system", which is either part of the documentation folder, included in the delivery of the device or available on request.
- ▷ If the interface specifications are missing or do not match with your specific DCU system, please contact your local representative of Sartorius Stedim Systems GmbH for support or write directly to

Sartorius Stedim Systems GmbH
Dept. FA (Product engineering, automation projects)
Robert-Bosch-Straße 5 – 7

D-34302 Guxhagen, Fed. Republic of Germany
phone +49.5665.407.0, fax +49.5665.407.2200

This page intentionally left blank.

Please provide these information to authorized personnel and the service only. If necessary, remove this page from the documentation and keep it on a safe place.

Certain system functions and settings which are limited to be done by qualified and authorized personnel only, are protected by the standard password system. This includes, for instance:

- ▷ the setting of controller parameters (e.g. PID) in the corresponding parameterization menus,
- ▷ in the main function "Settings"
 - the setting of process values "PV"
 - under "Manual Operation", e.g. for simulation, the setting of interface parameters of the digital and analog process inputs and outputs or of controllers.

In addition the submenu "Service" of the main menu "Settings" can only be accessed by a special service password [► fig. 8-1, -2], which will only be provided to authorized service staff.

At selection of a password-protected function or submenu, a corresponding key pad is displayed. The following passwords can be defined:

- ▷ Standard password, set in the factory: "[19]".
- ▷ Customized standard password,
as agreed on delivery: „[_____]" *).
- ▷ Service-Password: „[_____]" *)

*) delivered by extra mail or added to the [► Technical Documentation]

This page intentionally left blank.

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Sartorius Stedim Biotech GmbH
August-Spindler-Str. 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of
this publication may be reprinted
or translated in any form or by
any means without the prior
written permission of Sartorius
Stedim Biotech GmbH.
The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make changes
to the technology, features,
specifications and design of the
equipment without notice.

Status:
February 2014,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany

LevMixer® is a trademark or registered
trademark of ATMI, Inc. in the
United States, other countries or both.

Printed in the EU on paper bleached
without chlorine | W · KT
Publication No.: SKG6001-a140202
Ver. 02 | 2014